



Repose<sup>®</sup> Flex

Designed to reduce  
the risk of DVT<sup>1</sup>

repose<sup>®</sup>



Utilising immersion and envelopment, Repose Flex is designed to place the knee in slight flexion (5° to 10°) and allow complete heel offloading while effectively re-distributing the weight of the patient's leg across its surface<sup>2</sup>.

By placing the knee in slight flexion, the risk of Popliteal Vein Compression at the rear of the knee is reduced<sup>1</sup>.

Popliteal Vein Compression is a significant source of Deep Vein Thrombosis for patients undergoing surgery or prolonged bed rest<sup>3</sup>.



Suitable for:

- Patients identified as being at risk of developing a DVT
- Bedbound patients, or patients undergoing surgery
- Protecting patients at all levels of risk of developing all stages of pressure ulcer
- Treatment of pressure ulcers when used as part of a package of care
- Deflatable for convenient storage and transport
- Repose smart valve technology for optimum inflation levels
- Two years warranty as standard



Repose protects patients at all levels of risk



Repose is indicated for the treatment of all categories of pressure ulcers



Multi-Patient Use Product

Description	Size (L/W/H)	Order Code	NHS SC
<b>XS</b> Repose Flex & Pump - Extra Small	L 420mm x W 700mm x H 90mm	6750000	FER14150
<b>S</b> Repose Flex & Pump - Small	L 460mm x W 700mm x H 90mm	6750001	FER14151
<b>STD</b> Repose Flex & Pump - Standard	L 500mm x W 700mm x H 90mm	6750002	FER14152

## References

1. npuap-epuap-pppia pressure ulcer treatment & prevention 2014 quick reference guide
2. Data Report on File: Taylor D. Validation of product function, Aug 2018
3. Popliteal vein compression syndrome the MAIN cause of DVT, Dr David M Grosser MB BS FRACS FRCS (London), Jan 2017

For more information:

Call +44 (0)1495 235800

Email: [info@frontier-group.co.uk](mailto:info@frontier-group.co.uk)

Website: [www.frontier-group.co.uk](http://www.frontier-group.co.uk)



Frontier Medical Group, Block G Newbridge Road Industrial Estate, Blackwood, Caerphilly NP12 2YN United Kingdom.