



Repose Foot Protector Range



Repose Foot Protectors
Designed specifically to minimise
the risk of pressure damage to heels.

repose[®]

Repose Foot Protector

The Repose Foot Protector and Foot Protector Plus are designed specifically to minimise the risk of pressure damage to heels.

A significant proportion of pressure ulcers occur on the heels due to a combination of pressure and friction/shear and whilst support surfaces assist in reducing

Smart Valve Technology

The 'Smart valve' ensures Repose products can never be over-inflated.

Offloading

Repose Foot Protector and Foot Protector Plus are designed to completely off-load pressure from the heel.

A heel void area incorporated into the Repose Foot Protector provides complete heel offloading.



pressure, only a specific device such as Repose Foot Protector and Foot Protector Plus can ensure that pressure on the heel is completely offloaded and so minimising the risk of pressure damage.

Clinical studies have shown the Repose Foot Protectors to be effective in a wide variety of settings.



Repose products protect patients at all levels of risk

- ✔ VERY HIGH RISK
- ✔ HIGH RISK
- ✔ AT RISK

Repose Foot Protector Plus

The Repose Foot Protector Plus combines the qualities of Repose Foot Protector with the additional benefits of:

ID Label

Provides the option for product traceability plus also includes cleaning and washing instructions.

Ventilation & Visual Inspection Panel

The open heel void of Repose Foot Protector Plus increases ventilation of the foot and so contributes to patient comfort. The open heel void also allows for easier visual inspection of the heel for wound observation and monitoring.


Ventilation & Visual Inspection Panel



Magnaffix Adjustable Strap

Magnaffix is a unique, patented multi-adjustable air filled strap incorporating magnets. This retains the product quickly and easily but is also easy to clean and remove.

The magnets have been sealed within the air cell pocket barrier so that they do not come into contact with the patient's skin to minimise the risk of any pressure or skin tear damage.



Repose Foot Protector Extra Small

Product	1 x Foot Protector Extra Small (Pair) & 1 x Pump
Dimensions	L 290mm x W 200mm x H 180mm

New Size

Repose Foot Protector Standard

Product	1 x Foot Protector Standard (Pair) & 1 x Pump
Dimensions	L 380mm x W 200mm x H 230mm.

Repose Foot Protector Small

Product	1 x Foot Protector Small (Pair) and 1 x Pump
Dimensions	L 330mm x W 200mm x H 205mm.

New Size

Repose Foot Protector Large

Product	1 x Foot Protector Large (Pair) & 1 x Pump
Dimensions	L 510mm x W 190mm x H 330mm

Repose Foot Protector Plus Extra Small

Product	1 x Foot Protector Plus Extra Small (Pair) & 1 x Pump
Dimensions	L 290mm x W 200mm x H 180mm

New Size

Repose Foot Protector Plus Standard

Product	1 x Foot Protector Plus Standard (Pair) & 1 x Pump
Dimensions	L 380mm x W 200mm x H 230mm.

Repose Foot Protector Plus Small

Product	1 x Foot Protector Plus Small (Pair) and 1 x Pump
Dimensions	L 330mm x W 200mm x H 205mm


New Size

Sizing Guide

When selecting an appropriate size Repose Foot Protector, you should ensure to measure from the base of the heel to the mid-point of the back of the knee, whilst the patient is sitting.

This is shown in the diagram below. Correspond your measurements to the leg size chart to select an applicable product.

Leg Measurement	Product Size
460mm - 410mm	Extra Small
500mm - 450mm	Small
550mm - 490mm	Standard
690mm - 550mm	Long
Maximum calf circumference 330mm	



Order Codes

Description	Size (L/W/H)	Order Code	NHS SC
Repose Foot Protector Extra Small (Pair)	L 290mm x W 200mm x H 180mm	6531100	FER14103
Repose Foot Protector Small (Pair)	L 330mm x W 200mm x H 205mm	6521100	FER14105
Repose Foot Protector Standard (Pair)	L 380mm x W 200mm x H 230mm	6501100	FER731
Repose Foot Protector Long (Pair)	L 510mm x W 200mm x H 330mm	6021100	FER8147
Repose Foot Protector Plus Extra Small (Pair)	L 290mm x W 200mm x H 180mm	6551100	FER14107
Repose Foot Protector Plus Small (Pair)	L 330mm x W 200mm x H 205mm	6541100	FER14109
Repose Foot Protector Plus Standard (Pair)	L 380mm x W 200mm x H 230mm	6511100	ELY591

For more information:

Call +44 (0)1495 235800

Email: info@frontier-group.co.uk

Website: www.frontier-group.co.uk



Frontier Medical Group, Block G Newbridge Road Industrial Estate, Blackwood, Caerphilly NP12 2YN United Kingdom.

© Copyright Frontier Therapeutics Limited, 2018
All rights reserved. GP0041_A59_02

repose[®]