



Repose Mattress Overlay



Repose Mattress Overlay
Support surface pressure
redistribution overlay

repose[®]

Repose Mattress Overlay



Repose Double Mattress Overlay

Repose Mattress Overlay & Cushion



Repose products have a two year warranty



Increased weight limit 222kg*
*Repose Mattress Overlay Only



Durable

Increased weight limit of 222kg* makes Repose products even suitable for larger patients. The new two year warranty is supported by clinical evidence of prolonged use.¹



Smart

Smart valve technology provides reassurance that Repose products are always inflated to the optimum pressure and never over inflated.



Multi-patient Use

Repose can be cleaned and reused for multi-patient use, increasing the lifetime of the product.



Cost Effective

Multi-patient use and an increased warranty to two years, provides cost-effective pressure area care.



Comfortable

Silent and unobtrusive, Repose mattress overlays and cushions are breathable and with no hard edges, increasing patient comfort and preventing tissue damage from occurring.



Support

Dedicated sales and clinical support ensure you and your patients will be supported every step of the way.



Simple and effective

The Repose range adopts the most natural and effective method of pressure relief and redistribution by reducing the pressure exerted at the interface between the patient and the supporting surface.

This is achieved by immersion, maximising the contact surface area and redistributing pressure without compromising surrounding tissue.



Prevention and treatment

Repose cushion and mattress overlays are suitable for patients at all levels of risk including those at very high risk of developing a pressure injury.

- ✓ VERY HIGH RISK
- ✓ HIGH RISK
- ✓ AT RISK

Repose® is suitable for patients **AT RISK** to **VERY HIGH RISK**

Repose products can also be used for patients with existing pressure damage and providing treatment of grade 1 and 2 pressure ulcers.

Clinical efficacy

Repose is an in-expensive system with demonstrable clinical efficacy, supported with 20 years of clinical evidence.

The most recent study using the Repose mattress overlay, demonstrated that when used in combination with a high specification foam mattress there was a reduction in pressure ulcer incidence when compared to only using a high specification foam mattress (5.2% vs 22.2%).²



Simple
& Effective



Order Codes

Description	Size (L/W/H)	Order Code	NHS SC
Repose Mattress Overlay & Pump	L 1780mm x W 770mm	6101100	FER733
Repose Mattress Overlay, Cushion & Pump Set	Mattress Overlay: L 1780mm x W 770mm Cushion: L 450mm x W 450mm	6121100	FER732
Repose Double Mattress Overlay & Pump	L 1780mm x W 1400mm	6131100	-

Reference 1. Impact Of Durability On The Value Proposition Of A Low-Cost, Low-Tech Support Surface Product Range. Maria Pool, Lead Nurse Tissue Viability, Wolverhampton PCT & Ann Coghlan Wound Care Sister, South Staffordshire PCT

Reference 2. Pressure relief with visco-elastic foam or with combined static air overlay? A prospective, crossover randomised clinical trial in a Dutch Nursing HomeMartin van Leen, Wounds 2013; 25(10): 287-292

*222kg for a Repose mattress overlay and double mattress overlay only when used on a high specification foam mattress. 139kg weight limit for all Repose products plus the mattress overlay and double mattress overlay if used on a standard foam mattress. The Repose Babynest 3kg max.

For more information:

Call +44 (0)1495 235800

Email: info@frontier-group.co.uk

Website: www.frontier-group.co.uk



Frontier Medical Group, Block G Newbridge Road Industrial Estate,
Blackwood, Caerphilly NP12 2YN United Kingdom.

© Copyright Frontier Therapeutics Limited, 2018
All rights reserved. GP0076_A73_06

repose[®]