




Repose Wedge



Repose Wedge
Designed to completely
off-load heels

repose[™]

Repose Wedge



The Repose Wedge is designed to completely off-load heels for either pressure damage prevention or as part of the management of heel ulcers.

Repose Wedge is suitable for patients:

- Assessed as being up to very high risk of developing heel pressure damage and /or with up to Category 4 pressure ulcer
- Who are confused and may not be able to remove other off-loading devices when moving off the bed
- Who are restless and therefore may not tolerate other forms of heel protection or kick them off
- With wider calves or suffer from foot drop making regular foot protector boots unsuitable
- On their sides or tilted and where there may be the possibility of pressure damage from an inflation value



Repose Wedge is now available with a PU strap for easy cleaning and additional infection prevention.



Repose products protect patients at all levels of risk

- ✓ VERY HIGH RISK
- ✓ HIGH RISK
- ✓ AT RISK



Multi-Patient Use Product

| Description | Size (L/W/H) | Order Code | NPC |
|------------------------------|-------------------|------------|---------|
| Repose Wedge & Pump | L 690mm x W 450mm | 6711100 | FER8136 |
| Repose Wedge PU Strap & Pump | L 690mm x W 450mm | 6711111 | ELY619 |

For more information:

Call +44 (0)1495 235800

Email: info@frontier-group.co.uk

Website: www.frontier-group.co.uk



Frontier Medical Group, Block G Newbridge Road Industrial Estate, Blackwood, Caerphilly NP12 2YN United Kingdom.

© Copyright Frontier Therapeutics Limited, 2018
All rights reserved. GP0041_A41_08

