



Repose Product Range

Reactive Pressure Redistribution

repose[®]

Simple and effective



The Repose range adopts an effective method of pressure redistribution and offloading by reducing the pressure exerted at the interface between the patient and the supporting surface.

This is achieved by immersion, maximising the contact surface area and redistributing pressure without compromising surrounding tissue.

Benefits of Repose



Durable

A two year warranty is supported by clinical evidence of prolonged use (Poole, 2009)



Smart

Smart valve technology provides reassurance that Repose products are easily inflated to the optimum pressure and never over inflated



Multi-patient Use

Repose products are easy to clean in line with local policy and reusable for multi-patient use, increasing the lifetime of the product



Cost Effective

Multi-patient use and a two-year warranty, provides cost-effective pressure area care



Comfortable

Repose products are breathable and have no hard edges, increasing patient comfort and preventing tissue damage from occurring



Support

Dedicated sales and clinical support ensure clinicians and patients are supported every step of the care pathway



The Repose Range



Surface Solutions



Repose® Mattress Overlay
Product order code: 6101100
NHS SC: FER733
Dimensions: L 1780mm x W 770mm
Weight Limit: 222kg*



Repose® Double Mattress Overlay
Product order code: 6131100
Dimensions: L 1780mm x W 1350mm
Weight Limit: 222kg Per Patient*



Repose® Mattress Overlay & Cushion Set
Product order code: 6121100
NHS SC: FER732
Dimensions: Mattress L 1780mm x W 770mm and Cushion L 450mm x W 450mm
Weight Limit: Cushion - 139kg* Mattress - 222kg*



Repose® Trolley Mattress Overlay
Product order code: 6601100
NHS SC: FER8146
Dimensions: L 1780mm x W 550mm
Weight Limit: 139kg



Repose® Companion
Product order code: 6990003
NHS SC: FER8140
Dimensions: L 1970mm x W 880mm
Weight Limit: 139kg



Repose® Trolley Companion
Product order code: 6990050
Dimensions: L 2000mm x W 615mm
Weight Limit: 139kg



Seating Solutions



Repose® Cushion
Product order code: 6201100
NHS SC: FER734
Dimensions: L 450mm x W 450mm
Weight Limit: 139kg



Repose® Care-Sit & Cover 400mm
Product order code: 6350001
NHS SC: FER14054
Dimensions: Care-Sit L 950mm x W 400mm
 Care-Sit Cover L 950mm x W 400mm
Weight Limit: 139kg



Repose® Care-Sit & Cover 450mm
Product order code: 6351001
NHS SC: ELY589
Dimensions: Care-Sit L 950mm x W 450mm
 Care-Sit Cover L 950mm x W 450mm
Weight Limit: 139kg



Repose® Contur Overlay
Product order code: 6310001
NHS SC: FER12421
Dimensions: L 1760mm x W 520mm
Weight Limit: 139kg



Paediatric Solutions




Repose® Paediatric Mattress Overlay
Product order code: 6801100
NHS SC: FER8141
Dimensions: L 1290mm x W 670mm
Weight Limit: 139kg



Repose® Babynest
Product order code: 6801101
NHS SC: FER001
Dimensions: L 580 mm x W 330 mm x H 100 mm
Weight Limit: 3kg




Foot Care Solutions



New Sizes

Repose® Foot Protector Standard
Product order code: 6501100
NHS SC: FER731
Dimensions: L 380mm x W 200mm x H 230mm
Weight Limit: 139kg

Available sizes: XS S STD L



New Sizes

Repose® Foot Protector Plus Standard
Product order code: 6511100
NHS SC: ELY591
Dimensions: L 380mm x W 200mm x H 230mm
Weight Limit: 139kg

Available sizes: XS S STD



PU Strap

Repose® Wedge
Product order code: 6711100
NHS SC: FER8136
Dimensions: L 690mm x W 450mm
Weight Limit: 139kg

Repose® Wedge With PU Strap
Product order code: 6711111



Repose® Sole Protector
Product order code: 6930001
NHS SC: FER13070
Dimensions: L 860mm x H 370mm
Weight Limit: 139kg



Repose® Babytherm
Product order code: 6401100
NHS SC: FER8138
Dimensions: L 760mm x W 530mm
Weight Limit: 139kg

Accessories

+



Repose® Cushion Cover

Product order code: 6211100
NHS SC: TLC1020
Dimensions: L 520mm x W 520mm

Repose® Mattress Cover

Product order code: 6211110
Dimensions: L 1890mm x W 890mm



Repose® Refresh Cushion Cover

Product order code: 6231104
NHS SC: TLC1021
Dimensions: L450mm x W 450mm

Repose® Refresh Mattress Cover

Product order code: 6050101
NHS SC: FER8143
Dimensions: L 1780mm x W 770mm

Repose® Cushion Strap Extender

Product order code: 6000001
Dimensions: L 1000mm x W 25mm

References

1. Price. P. et al. 1999. Challenging the pressure sore paradigm. Journal of Wound Care, April, Vol 8, No 4.
2. Ballard Wilson, A. Pressure ulcer prevention audit: What are the benefits of doing it? Poster EPUAP Berlin, 2006.
3. Staffordshire County Council Loan Store 2007. Data on file.
4. Bale. S. Interim results of a randomised controlled trial of two pressure - relieving systems in the prevention of pressure sores in patients with fractured neck of femur. Presented at the 7th European Conference on Advances in Wound Management, Harrogate, 1997.
5. Osterbrink, J et al. Clinical evaluation of the effectiveness of a multimodal static pressure relieving device. 8th European Pressure Ulcer Advisory Panel Open Meeting, Aberdeen. May 5th – 7th, 2005.
6. NPUAP-EPUAP 2009 Pressure Ulcer Treatment QRG p19
7. Price P et al. 2003 The use of a new overlay mattress in patients with chronic pain: impact on sleep and self-reported pain. Clinical Rehabilitation 2003:17.
* Nimbus is a registered trademark of Huntleigh Technology PLC

*Mattress overlay and double mattress overlay when used on a high specification foam mattress 222kg.

*Mattress overlay and double mattress overlay when used on a non-high specification foam mattress 139kg.

For more information:

Call +44 (0)1495 235800

Email: info@frontier-group.co.uk

Website: www.frontier-group.co.uk



Frontier Medical Group, Block G Newbridge Road Industrial Estate,
Blackwood, Caerphilly NP12 2YN United Kingdom.

CE © Copyright Frontier Therapeutics Limited, 2018
All rights reserved. GP0041_A20_12

repose®